



## LEAK TEST APPARATUS

This apparatus is suitable for checking the sealing of Blister Packs and Strips. It is equipped with efficient powerful Vacuum Pump, which ensures the test in shortest possible time. A vacuum gauge is provided to indicate vacuum level. More details, as per enclosed illustrated leaflet. Unit is available in two models.

### Model: LTA-150.

Small Desiccator with 150 mm, diameter

### Model: LTA-300.

Large Desiccator with 300 mm, diameter



**Registered Office:** Survey No. 81, Plot No. 16-17, Junagadh Road, Near Lion's School, Jetpur - 360370, Gujarat, India.

**Branch Office:** 401, Vithal Bhavan, 1st Floor, Chira Bazaar J.S.S. Road, Mumbai - 400 002, India.

**Tel:** +91 022 66370811, +91 022 66370812 **Email:** info@alakhealthcare.com **Website:** www.alakhealthcare.com